

Cuscino in fibra cava siliconata Hollow fiber polyester cushion

pl 150N-150F

Codice Code	Misura (cm) Size (cm)
PL150/N	43X43X6
* PL150/F	43X43X6



* Versione con il foro
* Central hole version



Cover

7

Drill Jeans

Cuscino antidecubito in fibra cava siliconata, estremamente soffice ed elastico, grazie alla presenza di fibre contenute in una sacca di cotone sanforizzata che ne costituiscono l'imbottitura. Le stesse fibre sono rivestite da un sottile strato di silicone che consente un più facile scorrimento delle une sulle altre, evitando l'ammassamento delle stesse. Ciò garantisce una buona ripartizione delle pressioni di appoggio e un'ottima aerazione della cute dovuta anche alla trapuntatura. Disponibile anche nella versione con foro centrale, dotato di lacci per fissaggio a carrozzine.

Grado di rischio:
prevenzione.

Non sfoderabile

Anti-decubitus cushion manufactured in hollow fiber polyester. It is extremely fluffy and elastic, thanks to the filling fibers inside a sanforized cotton bag. The fibers, covered by a thin silicone layer, permit a better sliding between themselves, avoiding the clumping. This guarantees a good contact pressure distribution and an excellent skin ventilation (also due to the quilting). Available in a central-hole version, equipped with laces for the wheelchair fastening.

Risk level:
prevention.

The cover is not removable