

# pl 151

## Cuscino in fibra cava siliconata Hollow fiber polyester cushion

Codice Code	Misura (cm) Size (cm)
PL151	40X40X9



### Cover

**2** Bielastico  
Bistretch

Cuscino antidecubito in fibra cava siliconata, estremamente soffice ed elastico, grazie alla presenza di fibre contenute in una sacca di cotone sanforizzata che ne costituiscono l'imbottitura. Le stesse fibre sono rivestite da un sottile strato di silicone che consente un più facile scorrimento delle une sulle altre, evitando l'ammassamento delle stesse. Ciò garantisce una buona ripartizione delle pressioni di appoggio e un'ottima aerazione della cute dovuta anche alla trapuntatura. Dotato di cover sfoderabile per la sanificazione ed il lavaggio.

**Grado di rischio:**  
prevenzione.

Anti-decubitus cushion manufactured in hollow fiber polyester. It is extremely fluffy and elastic, thanks to the filling fibers inside a sanforized cotton bag.

The fibers, covered by a thin silicone layer, permit a better sliding between themselves, avoiding the clumping. This guarantees a good contact pressure distribution and an excellent skin ventilation (also due to the quilting). Equipped with removable cover for sanitizing and washing.

**Risk level:**  
prevention.

