


# RSFT/D

## Materasso in poliuretano espanso a densità differenziata Mattress in polyurethane foam with differentiated densities

Codice Code	Misura (cm) Size (cm)	Portanza funzionale (Kg) Max supported weight (Kg)
RSFT/D	195X85X17	160

*altre misure disponibili su richiesta  
other sizes available on request*



### Cover

- 1 Bielastico + antiscivolo (serigrafia su richiesta)  
Bistretch + antislip (serigraphy on request)

Materasso antidecubito preventivo in poliuretano espanso a densità differenziate, in cui lo strato superiore ad alta resilienza garantisce uno scarico ottimale delle forze di compressione sulla cute del paziente.

E' formato da una struttura portante con sponde laterali di contenimento e da un inserto centrale con canalizzazione orizzontale. Disponibile a sezione unica o a più sezioni e adattabile ai movimenti di tutti i letti articolati per garantire un alto comfort al paziente.

### Grado di rischio:

prevenzione e primo stadio.

Ventilated polyurethane foam anti-decubitus prevention mattress with differentiated densities where the high-resilience top layer guarantees an optimal sinking of the pressure forces of the patient skin.

Includes a bearing structure with lateral containment edges and a central polyurethane foam insert with horizontal root canal treatment. It is available both in singular-section and multiple-sections and is adaptable to every profiled beds' movement maintaining the high comfort for the patient.

### Risk level:

prevention and first level.