

## Cuscino anatomico con inserto in gel Anatomic cushion with gel inserts

pl 6000

Codice Code	Misura (cm) Size (cm)
PL6000/4040	40X40X8
PL6000/4242	42X42X8
PL6000/4545	45X45X8



### Cover

**2** Bielastico  
Bistretch

**5** Bielastico + 3D + antiscivolo  
Bistretch + 3D + antislip

Cuscino antidecubito costituito da una base in schiuma di poliuretano a "lenta memoria" su cui è presente un gel di poliuretano automodellante.

Le caratteristiche della base, di forma anatomica, consentono di sostenere il peso del paziente, adeguandosi alla sua morfologia ossea e quindi garantendo lo scarico ottimale delle pressioni di contatto. L'inserto in gel attenua ulteriormente le pressioni della zona ischiatica per un comfort migliore del paziente.

**Grado di rischio:**  
prevenzione.

Anti-decubitus cushion manufactured with a polyurethane memory foam base with an auto-shaping polyurethane gel insert on top.

The anatomic shape base characteristics allow to support the patient's weight, adapting to his osseous morphology therefore guaranteeing the optimal contact pressures sinking. The gel insert alleviates even more the sciatic zone pressures for the better patient comfort.

**Risk level:**  
prevention.