


Materasso in poliuretano espanso a lenta memoria Mattress in polyurethane memory foam

RLOW

Codice Code	Misura (cm) Size (cm)	Portanza funzionale (Kg) Max supported weight (Kg)
RLOW	195x85x14	160

*altre misure disponibili su richiesta
other sizes available on request*



Cover

- 1 Bielastico + antiscivolo (serigrafia inclusa)
Bistretch + antislip (serigraphy included)

Materasso antidecubito preventivo in poliuretano espanso a lenta memoria.

Formato da una struttura portante con sponde laterali di contenimento e da un inserto a lenta memoria nello strato superiore.

Adattabile ai movimenti di tutti i letti articolati. Disponibile liscio o con lavorazione ondulata, consente un'ottimale distribuzione delle pressioni garantendo uno scarico delle forze di compressione sulla cute, creando un comfort ideale per il paziente.

Grado di rischio:

prevenzione e primo stadio.

Polyurethane memory foam anti-decubitus prevention mattress.

Includes a bearing structure with lateral containment edges and a central memory foam insert. Adaptable to profiled beds' movement.

Available both in smooth and wavy shapes and guarantees an optimal sinking of the pressure forces on the patient skin creating an ideal comfort for the patient.

Risk level:

prevention and first level.