


TOPPER

Sovramaterasso in lenta memoria Overlay mattress in memory foam

Codice Code	Misura (cm) Size (cm)	Portanza funzionale (Kg) Max supported weight (Kg)
PLTOPPER	195X85X5	130

altre misure disponibili su richiesta
other sizes available on request



Cover

4 Bielastico
Bistretch

E' un sovr materasso antidecubito preventivo in poliuretano espanso a memoria di forma che garantisce un buon comfort del paziente.

Questa superficie può essere posta direttamente sopra al proprio materasso per migliorarne le caratteristiche. Grazie alla sua composizione che garantisce una migliore circolazione del sangue si ottiene un effetto ergonomico unico che ottimizza la postura e lo scarico della colonna vertebrale.

Grado di rischio:
prevenzione.

An overlay anti-decubitus prevention mattress made in memory foam which guarantees a nice comfort for the patient.

This surface can be placed directly over the mattress to improve its characteristics. Thanks to its composition a better blood flow is guaranteed, its ergonomic effect optimizes the spinal column's posture and the sinking.

Risk level:
prevention.

