

pl 5000

Cuscino in viscoelastico con inserto in gel Viscoelastic cushion with gel inserts

Codice Code	Misura (cm) Size (cm)
PL5000/4040	40X40X6
PL5000/4242	42X42X6
PL5000/4545	45X45X6

*altre misure disponibili su richiesta
other sizes available on request*



Cover

2 Bielastico
Bistretch

5 Bielastico + 3D + antiscivolo
Bistretch + 3D + antislip

Cuscino antidecubito in viscoelastico compatto costituito da uno strato di gel di poliuretano e uno strato di poliuretano espanso a lenta memoria. L'azione combinata dei due componenti garantisce una diminuzione delle pressioni di contatto e consente un significativo aumento della superficie di appoggio per evitare l'insorgenza di lesioni della cute del paziente con la superficie sottostante.

Grado di rischio:
prevenzione.

Anti-decubitus cushion manufactured with compact viscoelastic, composed by a polyurethane gel layer and an expanded memory foam polyurethane layer. Their combined action guarantees a decrease of the contact pressures and allows a big increase of the support surface therefore avoiding the onset of the patient skin lesions with the underlying surface.

Risk level:
prevention.

